


钛缆固定系统（无菌包装）

Titanium Cable Fixation System (Sterile Packaging)



系统特点

Character Of The System

- 由多根钛合金丝捻制而成，柔韧性好，具有出众的静态强度和更高的疲劳强度。
- 钛缆具有很强的抗应变能力，可以理想的拉紧和固定。
- 锁扣机制锁定，无需扭转线缆。
- 通过收紧器可精确控制张力。


- Made of multiple twisted titanium alloy wire, good flexibility, outstanding static strength and higher fatigue strength.
- Titanium cable has strong strain resistance capability, ideal for tensioning and fixing.
- Fasten buckle mechanism locking, no need of twisting wire cable.
- Control tension precisely with tightener.

适应症


Indications

髌骨骨折，尺骨鹰嘴骨折，肱骨大结节骨折，内踝骨折，肩锁关节脱位，股骨大转子骨折等骨性突起的一些撕脱骨折。

Patella fracture, olecranon fracture, humeral Large tuberosity fracture, medial condyle fracture, Acromioclavicular joint dislocation, femur trochanter fracture and other avulsion fractures of bony prominence.



鹰嘴骨折



髌骨骨折



人工关节假体翻修



趾骨骨折

钛缆固定系统（无菌包装）

Titanium Cable Fixation System (Sterile Packaging)

金属缆索 T771
Metal Titanium Cable



产品名称 Product Name	包装形式 Packaging Form	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
金属缆索 Metal Titanium Cable	无菌包装 Sterile Packaging	钛合金 Titanium Alloy	1.3	520	T7713520-S
			1.8	520	T7718520-S

锁紧扣 T772
Fasten Buckle



产品名称 Product Name	包装形式 Packaging Form	材质 Material	直径(mm) Diameter	产品编号 Product No.
锁紧扣 Fasten Buckle	无菌包装 Sterile Packaging	钛合金 Titanium Alloy	1.4	T7721408-S
			1.9	T7721908-S

张力带别针 T522
Tension Band Pins



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
张力带别针 Tension Band Pins	钛合金 Titanium Alloy	2.0	60	T5220060
			70	T5220070
			80	T5220080
			100	T5220100
			120	T5220120
		2.5	60	T5225060
			70	T5225070
			80	T5225080
			100	T5225100
			120	T5225120

无菌包装
Sterile Packaging



配套器械
Supporting Instruments

