

锁定钢板系列

Locking Plate Series

KUS 股骨近端外侧锁定钢板(左右) T4424
KUS-Proximal Lateral Femoral Locking Plate
(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
KUS 股骨近端外侧 锁定钢板(左右) KUS-Proximal Lateral Femoral Locking Plate (Left,Right)	纯钛 Pure Titanium	4	149	T4424004R/L
		5	169	T4424005R/L
		6	189	T4424006R/L
		8	229	T4424008R/L
		10	269	T4424010R/L
		12	309	T4424012R/L
备注：头部选用Φ6.5中空锁定螺钉，其余选用Φ5.0锁定螺钉和HA4.5螺钉(图为5孔) Note: used for Φ6.5 locking screw in head, other Φ5.0 locking screw and HA4.5 screw (Sample: 5 holes)				

KUH 股骨近端锁定钢板(左右) T4430
KUH-Proximal Femoral Locking Plate (Left,Right)




产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
KUH 股骨近端 锁定钢板(左右) KUH-Proximal Femoral Locking Plate (Left,Right)	纯钛 Pure Titanium	4	128	T4430004R/L
		5	148	T4430005R/L
		6	168	T4430006R/L
		8	208	T4430008R/L
		10	248	T4430010R/L
		12	288	T4430012R/L
备注：头部选用Φ6.5中空锁定螺钉，其余选用Φ5.0锁定螺钉和HA4.5螺钉(图为5孔) Note: used for Φ6.5 locking screw in head, other Φ5.0 locking screw and HA4.5 screw (Sample: 5 holes)				

KHS-Ⅰ 股骨近端锁定钢板(左右) T4425
KHS-Ⅰ Proximal Femoral Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
KHS-Ⅰ 股骨近端锁定钢板(左右) KHS-Ⅰ Proximal Femoral Locking Plate (Left,Right)	纯钛 Pure Titanium	4	T4425004R/L
备注：头部选用Φ7.0中空锁定螺钉，其余选用Φ5.0锁定螺钉 Note: used for Φ7.0 locking screw in head, other Φ5.0 locking screw			

KHS-Ⅱ 股骨近端加压锁定板(左右) T4426/T4427
KHS-Ⅱ Proximal Femoral Locking Compression Plate
(Left,Right)



视频在线播放

产品名称 Product Name	材质 Material	规格 Specification	产品编号 Product No.
KHS-Ⅱ 股骨近端加压锁定板(左右) KHS-Ⅱ Proximal Femoral Locking Compression Plate (Left,Right)	纯钛 Pure Titanium	大号3孔 Large 3 Holes	T4426003R/L
		小号3孔 Small 3 Holes	T4427003R/L

KHS-Ⅱ 股骨颈固定螺钉 T4907
KHS-Ⅱ Femoral Neck Fixing Screw



产品名称 Product Name	材质 Material	螺纹直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
KHS-Ⅱ 股骨颈固定螺钉 KHS-Ⅱ Femoral Neck Fixing Screw	钛合金 Titanium Alloy	7.0	60-115 (每5mm为一档) (The specification is 5mm)	T49070060 -T49070115

KHS-Ⅱ 支头螺丝 T4910
KHS-Ⅱ Nut



产品名称 Product Name	材质 Material	螺纹直径(mm) Diameter	产品编号 Product No.
KHS-Ⅱ 支头螺丝 KHS-Ⅱ Nut	钛合金 Titanium Alloy	9.0	T4910000

锁定钢板系列

Locking Plate Series

130°DHS 鹅头钉锁定钢板 T4500
130°DHS Compression Hip Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
130° DHS 鹅头钉 锁定钢板 130° DHS Compression Hip Locking Plate	纯钛 Pure Titanium	4	125	T4500004
		5	146	T4500005
		6	165	T4500006
		8	203	T4500008
备注: 选用Φ5.0锁定螺钉, HA4.5螺钉, HB6.5螺钉 (图例为5孔) Note: used for Φ5.0 locking screw, HA4.5 screw, HB6.5 screw (Sample: 5 holes)				

股骨远端外侧锁定钢板(左右) T4411
Distal Lateral Femoral Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
股骨远端外侧 锁定钢板(左右) Distal Lateral Femoral Locking Plate (Left,Right)	纯钛 Pure Titanium	5	166	T4411005R/L
		7	206	T4411007R/L
		9	246	T4411009R/L
		11	285	T4411011R/L
		13	325	T4411013R/L
备注: 选用Φ5.0锁定螺钉, HA4.5螺钉 (图例为9孔) Note: used for Φ5.0 locking screw, HA4.5 screw (Sample: 9 holes)				

胫骨近端外侧锁定钢板Ⅳ型(左右) T4700
Proximal Lateral Tibial Locking Plate-Ⅳ (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
胫骨近端外侧锁定 钢板Ⅳ型(左右) Proximal Lateral Tibial Locking Plate-Ⅳ (Left,Right)	纯钛 Pure Titanium	4	113	T4700004R/L
		5	129	T4700005R/L
		7	161	T4700007R/L
		9	193	T4700009R/L
		11	225	T4700011R/L
		13	257	T4700013R/L
备注：选用Φ5.0锁定螺钉，HA4.5螺钉（图例为7孔） Note: used for Φ5.0 locking screw, HA4.5 screw (Sample: 7 holes)				

KSL 胫骨近端外侧锁定钢板(左右) T4711
KSL Proximal Lateral Tibial Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
KSL 胫骨近端外侧 锁定钢板(左右) KSL Proximal Lateral Tibial Locking Plate (Left,Right)	纯钛 Pure Titanium	5	87	T4711005R/L
		7	116	T4711007R/L
		9	144	T4711009R/L
		11	172	T4711011R/L
		13	200	T4711013R/L
		15	228	T4711015R/L
备注：选用Φ3.5锁定螺钉，HA3.5螺钉（图例为9孔） Note: used for Φ3.5 locking screw, HA3.5 screw (Sample: 9 holes)				

3.5胫骨近端内侧锁定钢板(左右) T4770
3.5 Proximal Medial Tibial Locking Plate
(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
3.5胫骨近端内 侧锁定钢板 (左右) 3.5 Proximal Medial Tibial Locking Plate (Left,Right)	纯钛 Pure Titanium	5	101	T4770005R/L
		7	129	T4770007R/L
		9	157	T4770009R/L
		11	185	T4770011R/L
备注：选用Φ3.5锁定螺钉(图例为7孔) Note: used for Φ3.5 locking screw(Sample: 7 holes)				

下肢T型重建锁定钢板(左右) T4730
T Reconstruction Locking Plate (Lower Limb)
(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
下肢T型重建 锁定钢板(左右) T Reconstruction Locking Plate (Lower Limb) (Left,Right)	纯钛 Pure Titanium	7	91	T4730007R/L
		9	123	T4730009R/L
		11	155	T4730011R/L
备注: 选用Φ5.0锁定螺钉, HA4.5螺钉 (图例为9孔) Note: used for Φ5.0 locking screw, HA4.5 screw (Sample: 9 holes)				

锁定钢板系列

Locking Plate Series

下肢L型重建锁定钢板(左右) T4740
L Reconstruction Locking Plate (Lower Limb)
(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
下肢L型重建 锁定钢板(左右) L Reconstruction Locking Plate (Lower Limb) (Left,Right)	纯钛 Pure Titanium	7	101	T4740007R/L
		9	133	T4740009R/L
		11	165	T4740011R/L
备注: 选用Φ5.0锁定螺钉, HA4.5螺钉 (图例为7孔) Note: used for Φ5.0 locking screw, HA4.5 screw (Sample: 7 holes)				

胫骨平台后柱锁定钢板(左右) T4750
Tibial Posterior Column Platform Locking
Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
胫骨平台后柱 锁定钢板(左右) Tibial Posterior Column Platform Locking Plate (Left,Right)	纯钛 Pure Titanium	3	65	T4750003R/L
		5	91	T4750005R/L
		7	117	T4750007R/L
备注：选用Φ3.5锁定螺钉（图例为7孔） Note: used for Φ3.5 locking screw (Sample: 7 holes)				

胫骨平台内柱锁定钢板(左右) T4760
Tibial Interior Column Platform Locking Plate
(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
胫骨平台内柱 锁定钢板(左右) Tibial Interior Column Platform Locking Plate (Left,Right)	纯钛 Pure Titanium	4	61	T4760004R/L
		6	87	T4760006R/L
		8	113	T4760008R/L
备注：选用Φ3.5锁定螺钉（图例为8孔） Note: used for Φ3.5 locking screw (Sample: 8 holes)				

胫骨远端前外侧锁定钢板(左右) T4630
Distal Anterolateral Tibial Locking Plate
(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
胫骨远端前外侧 锁定钢板(左右) Distal Anterolateral Tibial Locking Plate (Left,Right)	纯钛 Pure Titanium	5	104	T4630005R/L
		7	130	T4630007R/L
		9	156	T4630009R/L
		11	182	T4630011R/L
		13	208	T4630013R/L
备注：选用Φ3.5锁定螺钉，HA3.5螺钉（图例为11孔） Note: used for Φ3.5 locking screw, HA3.5 screw (Sample: 11 holes)				

胫骨远端外侧重建锁定钢板(左右) T4600
Distal Lateral Tibial Reconstruction Locking Plate
(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
胫骨远端外侧重建 锁定钢板(左右) Distal Lateral Tibial Reconstruction Locking Plate (Left,Right)	纯钛 Pure Titanium	5	125	T4600005R/L
		7	161	T4600007R/L
		9	197	T4600009R/L
		11	233	T4600011R/L
		13	269	T4600013R/L
备注：选用Φ5.0锁定螺钉（图例为7孔） Note: used for Φ5.0 locking screw (Sample: 7 holes)				

锁定钢板系列

Locking Plate Series

胫骨远端内侧锁定钢板(左右) T4620
Distal Medial Tibial Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
胫骨远端内侧 锁定钢板(左右) Distal Medial Tibial Locking Plate (Left,Right)	纯钛 Pure Titanium	5	139	T4620005R/L
		7	167	T4620007R/L
		9	195	T4620009R/L
		11	223	T4620011R/L

备注：选用Φ3.5锁定螺钉，HA3.5螺钉（图例为7孔）
Note: used for Φ3.5 locking screw, HA3.5 screw (Sample: 7 holes)

骨盆重建锁定钢板 T4790
Pelvic Reconstruction Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
骨盆重建锁定钢板 Pelvic Reconstruction Locking Plate	纯钛 Pure Titanium	6	T4790006
		8	T4790008
		10	T4790010
		12	T4790012
		14	T4790014

备注：选用Φ3.5锁定螺钉（图例为12孔）
Note: used for Φ3.5 locking screw (Sample: 12 holes)

2.7 腓骨远端重建锁定钢板(左右) T4405
2.7 Distal Fibula Reconstruction Locking Plate
(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
2.7 腓骨远端重建 锁定钢板(左右) 2.7 Distal Fibula Reconstruction Locking Plate (Left,Right)	纯钛 Pure Titanium	4	68	T4405004R/L
		6	88	T4405006R/L
		8	108	T4405008R/L
		10	128	T4405010R/L
		12	148	T4405012R/L

备注：选用Φ2.7锁定螺钉，HA2.7螺钉（图例为6孔）
Note: used for Φ2.7 locking screw and HA2.7 screw (Sample: 6 holes)

肱骨近端锁定钢板Ⅲ型(左右) T4230
Proximal Humerus Locking Plate-Ⅲ (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
肱骨近端锁定 钢板Ⅲ型(左右) Proximal Humerus Locking Plate-Ⅲ (Left,Right)	纯钛 Pure Titanium	4	103	T4230004R/L
		5	117	T4230005R/L
		6	131	T4230006R/L
		8	159	T4230008R/L
		10	187	T4230010R/L
		12	215	T4230012R/L

备注：选用Φ3.5锁定螺钉，HA3.5螺钉，HB4.0螺钉（图例为6孔）
Note: used for Φ3.5 locking screw, HA3.5 screw, HB4.0 screw (Sample: 6 holes)

肱骨近端锁定钢板Ⅱ型 T4220
Proximal Humerus Locking Plate-Ⅱ



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
肱骨近端锁定 钢板Ⅱ型 Proximal Humerus Locking Plate-Ⅱ	纯钛 Pure Titanium	3	96	T4220003
		4	108	T4220004
		5	120	T4220005
		6	132	T4220006
		8	156	T4220008
		10	180	T4220010
		12	204	T4220012

备注：选用Φ3.5锁定螺钉，HA3.5螺钉，HB4.0螺钉（图例为6孔）
Note: used for Φ3.5 locking screw, HA3.5 screw, HB4.0 screw (Sample: 6 holes)

锁定钢板系列

Locking Plate Series

肱骨远端外侧锁定钢板(左右) T4270

Distal Lateral Humerus Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
肱骨远端外侧 锁定钢板(左右) Distal Lateral Humerus Locking Plate (Left,Right)	纯钛 Pure Titanium	5	100	T4270005R/L
		7	124	T4270007R/L
		9	148	T4270009R/L
		11	172	T4270011R/L
备注：头部选用Φ2.7锁定螺钉，其余选用Φ3.5锁定螺钉和HA3.5螺钉（图例为7孔） Note: used for Φ2.7 locking screw in head, other Φ3.5 locking screw and HA3.5 screw (Sample: 7 holes)				

肱骨远端内侧锁定钢板(左右) T4280

Distal Medial Humerus Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
肱骨远端内侧 锁定钢板(左右) Distal Medial Humerus Locking Plate (Left,Right)	纯钛 Pure Titanium	4	82	T4280004R/L
		5	94	T4280005R/L
		7	118	T4280007R/L
		9	142	T4280009R/L
备注：头部选用Φ2.7锁定螺钉，其余选用Φ3.5锁定螺钉和HA3.5螺钉(图例为7孔) Note: used for Φ2.7 locking screw in head, other Φ3.5 locking screw and HA3.5 screw (Sample: 7 holes)				

锁骨远端锁定钢板(左右) T4381

Distal Clavicle Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
锁骨远端锁定 钢板(左右) Distal Clavicle Locking Plate (Left,Right)	纯钛 Pure Titanium	4	72	T4381004R/L
		5	84	T4381005R/L
		6	96	T4381006R/L
		7	108	T4381007R/L
		8	119	T4381008R/L
		10	141	T4381010R/L
备注：头部选用Φ2.7锁定螺钉，其余选用Φ3.5锁定螺钉（图例为4孔） Note: used for Φ2.7 locking screw in head, other Φ3.5 locking screw (Sample: 4 holes)				

3.0 锁骨S型锁定钢板(左右) T4382

3.0 S Clavicle Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
3.0 锁骨S型 锁定钢板(左右) 3.0 S Clavicle Locking Plate (Left,Right)	纯钛 Pure Titanium	6	70	T4382006
		7	82	T4382007R/L
		8	94	T4382008R/L
		9	105	T4382009R/L
		10	116	T4382010R/L
		12	137	T4382012R/L
备注：选用Φ3.5锁定螺钉（图例为7孔） Note: used for Φ3.5 locking screw (Sample: 7 holes)				

锁定经皮锁骨钢板(左右) T4383

Percutaneous Clavicle Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
锁定经皮锁骨 钢板(左右) Percutaneous Clavicle Locking Plate (Left,Right)	纯钛 Pure Titanium	7	82	T4383007R/L
		8	94	T4383008R/L
		9	105	T4383009R/L
		10	116	T4383010R/L
		12	137	T4383012R/L
备注：选用Φ3.5锁定螺钉（图例为9孔） Note: used for Φ3.5 locking screw (Sample: 9 holes)				

锁定钢板系列

Locking Plate Series

肩锁钩锁定钢板(左右) T4390、T4391
Clavicle Hook Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	高度(mm) Height	孔数 Holes	长度(mm) Length	产品编号 Product No.
肩锁钩锁定 钢板(左右) Clavicle Hook Locking Plate (Left,Right)	纯钛 Pure Titanium	14	6	57	T4390006R/L
			7	69	T4390007R/L
			8	81	T4390008R/L
			9	93	T4390009R/L
			10	105	T4390010R/L
		17	6	57	T4391006R/L
			7	69	T4391007R/L
			8	81	T4391008R/L
			9	93	T4391009R/L
			10	105	T4391010R/L

备注：选用Φ3.5锁定螺钉（图例为7孔）
Note: used for Φ3.5 locking screw (Sample: 7 holes)

鹰嘴钩重建锁定钢板(左右) T4310
Olecranon Hook Reconstruction Locking
Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
鹰嘴钩重建 锁定钢板(左右) Olecranon Hook Reconstruction Locking Plate (Left,Right)	纯钛 Pure Titanium	4	87	T4310004R/L
		6	113	T4310006R/L
		8	139	T4310008R/L

备注：选用Φ3.5锁定螺钉（图例为5孔）
Note: used for Φ3.5 locking screw (Sample: 5 holes)

尺骨近端锁定钢板(左右) T4370
Proximal Ulna Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
尺骨近端锁 定钢板(左右) Proximal Ulna Locking Plate (Left,Right)	纯钛 Pure Titanium	4	119	T4370004R/L
		6	145	T4370006R/L
		8	171	T4370008R/L

备注：选用Φ3.5锁定螺钉（图例为5孔）
Note: used for Φ3.5 locking screw (Sample: 5 holes)

2.7 桡骨头锁定钢板(左右) T4322
2.7 Radius Head Locking Plate
(Left,Right)



视频在线播放

产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
2.7 桡骨头锁定钢板(左右) 2.7 Radius Head Locking Plate(Left,Right)	纯钛 Pure Titanium	3	40	T4322003R/L
		4	48	T4322004R/L
		5	56	T4322005R/L

备注：选用Φ2.7锁定螺钉（图例为5孔）
Note: used for Φ2.7 locking screw (Sample: 5 holes)

锁定钢板系列

Locking Plate Series

小T型重建锁定钢板 T4330
Mini-T Reconstruction Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
小T型重建锁定钢板 Mini-T Reconstruction Locking Plate	纯钛 Pure Titanium	6	43	T4330006
		7	54	T4330007
		8	65	T4330008
		9	76	T4330009

备注: 选用Φ3.5锁定螺钉 (图例为7孔)
Note: used for Φ3.5 locking screw (Sample: 7 holes)

桡骨斜T型重建锁定钢板(左右) T4340
Oblique T Radius Reconstruction Locking Plate
(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
桡骨斜T型重建 锁定钢板(左右) Oblique T Radius Reconstruction Locking Plate (Left,Right)	纯钛 Pure Titanium	6	45	T4340006R/L
		7	56	T4340007R/L
		8	67	T4340008R/L
		9	78	T4340009R/L

备注: 选用Φ3.5锁定螺钉 (图例为7孔)
Note: used for Φ3.5 locking screw (Sample: 7 holes)

2.7 桡骨远端支持锁定钢板(大号)(左右) T4320
2.7 Distal Radius Buttress Locking Plate
(Large)(Left,Right)



视频在线播放

产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
2.7 桡骨远端支持 锁定钢板(大号)(左右) 2.7 Distal Radius Buttress Locking Plate (Large)(Left,Right)	纯钛 Pure Titanium	4	55	T4320004R/L
		5	63	T4320005R/L
		7	79	T4320007R/L
		9	95	T4320009R/L

备注: 选用Φ2.7锁定螺钉, HA2.7螺钉 (图例为5孔)
Note: used for Φ2.7 locking screw, HA2.7 screw (Sample: 5 holes)

2.7 桡骨远端支持锁定钢板(小号)(左右) T4324
2.7 Distal Radius Buttress Locking Plate
(Mini)(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
2.7 桡骨远端支持 锁定钢板(小号)(左右) 2.7 Distal Radius Buttress Locking Plate (Mini)(Left,Right)	纯钛 Pure Titanium	4	55	T4324004R/L
		5	63	T4324005R/L
		7	79	T4324007R/L
		9	95	T4324009R/L

备注: 选用Φ2.7锁定螺钉, HA2.7螺钉 (图例为4孔)
Note: used for Φ2.7 locking screw, HA2.7 screw (Sample: 4 holes)

2.7 桡骨远端背侧锁定钢板(左右) T4321
2.7 Distal Back Radius Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.7 桡骨远端背侧 锁定钢板(左右) 2.7 Distal Back Radius Locking Plate (Left,Right)	纯钛 Pure Titanium	14	T4321014R/L

备注: 选用Φ2.7锁定螺钉 (图例为14孔)
Note: used for Φ2.7 locking screw (Sample: 14 holes)

锁定钢板系列

Locking Plate Series

跟骨锁定钢板 I 型(左右) T4400
Calcaneus Locking Plate- I (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
跟骨锁定钢板 I 型(左右) Calcaneus Locking Plate- I (Left,Right)	纯钛 Pure Titanium	8	49	T4400008R/L
		10	59	T4400010R/L
备注：选用Φ3.5锁定螺钉（图例为8孔） Note: used for Φ3.5 locking screw (Sample: 8 holes)				

跟骨锁定钢板 II 型(大号)(左右) T4402
Calcaneus Locking Plate- II (Max)(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
跟骨锁定钢板 II 型 (大号)(左右) Calcaneus Locking Plate- II (Max) (Left,Right)	纯钛 Pure Titanium	12	67	T4402012R/L
备注: 选用Φ3.5锁定螺钉 Note: used for Φ3.5 locking screw				

跟骨锁定钢板 II 型(中号)(左右) T4401
Calcaneus Locking Plate- II (Mid)(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
跟骨锁定钢板 II 型 (中号)(左右) Calcaneus Locking Plate- II (Mid) (Left,Right)	纯钛 Pure Titanium	12	58	T4401012R/L
备注: 选用Φ3.5锁定螺钉 Note: used for Φ3.5 locking screw				

跟骨锁定钢板 II 型(小号)(左右) T4401
Calcaneus Locking Plate- II (Mini)(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
跟骨锁定钢板 II 型 (小号)(左右) Calcaneus Locking Plate- II (Mini) (Left,Right)	纯钛 Pure Titanium	8	51	T4401008R/L
备注: 选用Φ3.5锁定螺钉 Note: used for Φ3.5 locking screw				

跟骨锁定钢板 II 型(短尾)(左右) T4404
Calcaneus Locking Plate- II (short tail)(Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
跟骨锁定钢板Ⅱ型(短尾)(左右) Calcaneus Locking Plate-Ⅱ (short tail) (Left,Right)	纯钛 Pure Titanium	11	59	T4404011R/L
		13	69	T4404013R/L
备注：选用Φ3.5锁定螺钉，HB4.0螺钉（图例为13孔） Note: used for Φ3.5 locking screw, HB4.0 screw (Sample: 13 holes)				

锁定钢板系列

Locking Plate Series

重建锁定钢板 T4040
Reconstruction Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
重建锁定钢板 Reconstruction Locking Plate	纯钛 Pure Titanium	5	63	T4040005
		6	76	T4040006
		7	89	T4040007
		8	102	T4040008
		9	115	T4040009
		10	128	T4040010
		12	154	T4040012
备注：选用Φ3.5锁定螺钉（图例为7孔） Note: used for Φ3.5 locking screw (Sample: 7 holes)				

1/3 锁定钢板 T4100
One-third Tubular Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
1/3 锁定钢板 One-third Tubular Locking Plate	纯钛 Pure Titanium	5	64	T4100005
		6	80	T4100006
		7	92	T4100007
		8	104	T4100008
		10	128	T4100010
备注: 选用Φ3.5锁定螺钉, HA3.5螺钉 (图例为6孔) Note: used for Φ3.5 locking screw, HA3.5 screw (Sample: 6 holes)				

尺桡骨锁定钢板 T4090
Ulna Radius Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
尺桡骨锁定钢板 Ulna Radius Locking Plate	纯钛 Pure Titanium	5	69	T4090005
		6	85	T4090006
		7	91	T4090007
		8	109	T4090008
		10	133	T4090010
		12	157	T4090012
备注：选用Φ3.5锁定螺钉，HA3.5螺钉（图例为7孔） Note: used for Φ3.5 locking screw, HA3.5 screw (Sample: 7 holes)				

肱骨锁定钢板 T4050
Humeral Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
肱骨锁定钢板 Humeral Locking Plate	纯钛 Pure Titanium	6	105	T4050006
		7	119	T4050007
		8	133	T4050008
		9	147	T4050009
		10	161	T4050010
		12	189	T4050012
		14	217	T4050014
备注：选用Φ3.5锁定螺钉，HA3.5螺钉（图例为7孔） Note: used for Φ3.5 locking screw, HA3.5 screw (Sample: 7 holes)				

胫骨锁定钢板 T4071
Tibial Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
胫骨锁定钢板 Tibial Locking Plate	纯钛 Pure Titanium	6	127	T4071006
		7	143	T4071007
		8	159	T4071008
		10	191	T4071010
		12	223	T4071012
		14	255	T4071014
备注: 选用Φ5.0锁定螺钉, HA4.5螺钉 (图例为7孔) Note: used for Φ5.0 locking screw, HA4.5 screw (Sample: 7 holes)				

股骨锁定钢板 T4083
Femoral Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	长度(mm) Length	产品编号 Product No.
股骨锁定钢板 Femoral Locking Plate	纯钛 Pure Titanium	8	195	T4083008
		9	215	T4083009
		10	235	T4083010
		12	275	T4083012
		14	315	T4083014
		16	355	T4083016
备注：选用Φ5.0锁定螺钉，HA4.5螺钉（图例为10孔） Note: used for Φ5.0 locking screw, HA4.5 screw (Sample: 10 holes)				

锁定钢板 2.0 系列

2.0 Locking Plate Series

2.0 T型锁定钢板 头部2孔 T41201

2.0 T-Type Locking Plate (Two Holes In The Head)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.0 T型锁定钢板 头部2孔 2.0 T-Type Locking Plate (Two Holes In The Head)	纯钛 Pure Titanium	5/6/7/8/9	T4120105-T4120109
备注：选用Φ2.0锁定螺钉(三角型)，HA2.0螺钉(三角型) (图例为7孔) Note: used for Φ2.0 locking screw (Triangle), HA2.0 screw (Triangle) (Sample: 7 holes)			

2.0 T型锁定钢板 头部3孔 T41202

2.0 T-Type Locking Plate (Three Holes In The Head)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.0 T型锁定钢板 头部3孔 2.0 T-Type Locking Plate (Three Holes In The Head)	纯钛 Pure Titanium	6/7/8/9	T4120206-T4120209
备注：选用Φ2.0锁定螺钉(三角型)，HA2.0螺钉(三角型) (图例为8孔) Note: used for Φ2.0 locking screw (Triangle), HA2.0 screw (Triangle) (Sample: 8 holes)			

2.0 踝部锁定钢板 T41204

2.0 Ankle Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.0 踝部锁定钢板 2.0 Ankle Locking Plate	纯钛 Pure Titanium	5/6/7/8/9	T4120405-T4120409
备注：选用Φ2.0锁定螺钉(三角型)，HA2.0螺钉(三角型) (图例为9孔) Note: used for Φ2.0 locking screw (Triangle), HA2.0 screw (Triangle) (Sample: 9 holes)			

2.0 L型锁定钢板(左右) T41206

2.0 L-Type Locking Plate (Left,Right)



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.0 L型锁定钢板(左右) 2.0 L-Type Locking Plate (Left,Right)	纯钛 Pure Titanium	5/6/7/8	T4120605-T4120608R/L
备注：选用Φ2.0锁定螺钉(三角型)，HA2.0螺钉(三角型) (图例为8孔) Note: used for Φ2.0 locking screw (Triangle), HA2.0 screw (Triangle) (Sample: 8 holes)			

2.0 直型锁定钢板 T41200

2.0 Straight Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.0 直型锁定钢板 2.0 Straight Locking Plate	纯钛 Pure Titanium	4/5/6/7/8/9/10	T4120004-T4120010
备注：选用Φ2.0锁定螺钉(三角型) (图例为6孔) Note: used for Φ2.0 locking screw (Triangle) (Sample: 6 holes)			

三角形锁定螺钉 T480

Triangle Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
三角形锁定螺钉 Triangle Locking Screw	钛合金 Titanium Alloy	2.0	6、8、10、12、 14、16、18、20、 22	T4802006- T4802022

金属接骨螺钉(三角型) T541

Metal Bone Fracture Screw (Triangle)



钉帽形式 Head of Nail Type	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
三角型 Triangle	钛合金 Titanium Alloy	2.0	6、8、10、12、 14、16、18、20、 22	T5412006- T5412022

锁定钢板 2.5 系列

2.5 Locking Plate Series

2.5 T型锁定钢板 头部2孔 T41251

2.5 T-Type Locking Plate (Two Holes In The Head)



2.5 T型锁定钢板 头部3孔 T41252

2.5 T-Type Locking Plate (Three Holes In The Head)



2.5 踝部锁定钢板 T41254

2.5 Ankle Locking Plate



2.5 L型锁定钢板(左右) T41256

2.5 L-Type Locking Plate (Left,Right)



2.5 直型锁定钢板 T41250

2.5 Straight Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.5 T型锁定钢板 头部2孔 2.5 T-Type Locking Plate (Two Holes In The Head)	纯钛 Pure Titanium	5/6/7/8/9	T4125105-T4125109
备注: 选用Φ2.5锁定螺钉, HA2.5螺钉 (图例为7孔) Note: used for Φ2.5 locking screw, HA2.5 screw (Sample: 7 holes)			

产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.5 T型锁定钢板 头部3孔 2.5 T-Type Locking Plate (Three Holes In The Head)	纯钛 Pure Titanium	6/7/8/9/10	T4125206-T4125210
备注: 选用Φ2.5锁定螺钉, HA2.5螺钉 (图例为8孔) Note: used for Φ2.5 locking screw, HA2.5 screw (Sample: 8 holes)			

产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.5 踝部锁定钢板 2.5 Ankle Locking Plate	纯钛 Pure Titanium	5/6/7/8/9	T4125405-T4125409
备注: 选用Φ2.5锁定螺钉, HA2.5螺钉 (图例为9孔) Note: used for Φ2.5 locking screw, HA2.5 screw (Sample: 9 holes)			

产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.5 L型锁定钢板(左右) 2.5 L-Type Locking Plate (Left,Right)	纯钛 Pure Titanium	5/6/7/8/9	T4125605-T4125609R/L
备注: 选用Φ2.5锁定螺钉, HA2.5螺钉 (图例为8孔) Note: used for Φ2.5 locking screw, HA2.5 screw (Sample: 8 holes)			


产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.5 直型锁定钢板 2.5 Straight Locking Plate	纯钛 Pure Titanium	4/5/6/7/8/10	T4125004-T4125010
备注: 选用Φ2.5锁定螺钉 (图例为6孔) Note: used for Φ2.5 locking screw (Sample: 6 holes)			

锁定钢板 2.7 系列

2.7 Locking Plate Series

2.7 H型锁定钢板 T41277

2.7 H-Type Locking Plate



2.7 X型锁定钢板 T41278

2.7 X-Type Locking Plate



产品名称 Product Name	材质 Material	孔数 Holes	产品编号 Product No.
2.7 H型锁定钢板 2.7 H-Type Locking Plate	纯钛 Pure Titanium	5	T4127705
		8	T4127708
备注：选用Φ2.7锁定螺钉 Note: used for Φ2.7 locking screw			

产品名称 Product Name	材质 Material	宽度(mm) Width	长度(mm) Length	产品编号 Product No.
2.7 X型锁定钢板 2.7 X-Type Locking Plate	纯钛 Pure Titanium	8	16	T4127816
		12	18	T4127818
		14	24	T4127824
			30	T4127830
备注：选用Φ2.7锁定螺钉 Note: used for Φ2.7 locking screw				

锁定螺钉系列

Locking Screw Series

三角型锁定螺钉 T480 Triangle Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
三角型锁定螺钉 Triangle Locking Screw	钛合金 Titanium Alloy	2.0	6、8、10、12、 14、16、18、20、 22	T4802006- T4802022

梅花型锁定螺钉 T481/T482 Torx Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型锁定螺钉 Torx Locking Screw	钛合金 Titanium Alloy	2.5	6、8、10、12、 14、16、18、20、 22、24、26、28、 30	T4812506- T4812530
		2.7	6、8、10、12、 14、16、18、20、 22、24、26、28、	T4822706- T4822728
			30、35、40、45、 50、55、60	T4822730- T4822760

梅花型锁定螺钉 T483/T484 Torx Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型锁定螺钉 Torx Locking Screw	钛合金 Titanium Alloy	3.5	10、12、14、16、 18、20、22、24、 26、28、30、32、 34、36、38、40、 42、44、46、48、	T4833510- T4833548
			50、55、60、65、 70、75、80	T4833550- T4833580
		5.0	14、16、18、20、 22、24、26、28、 30、32、34、36、 38、40、42、44、 46、48、	T4845016- T4845048
			50、55、60、65、 70、75、80、85、 90、95、100	T4845050- T4845100

梅花型中空锁定螺钉 T487 Torx Hollow Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型中空 锁定螺钉 Torx Hollow Locking Screw	钛合金 Titanium Alloy	6.5	30、35、40、45、 50、55、60、65、 70、75、80、85、 90、95、100、 105、110	T48765030- T48765110

梅花型中空锁定螺钉 T488 Torx Hollow Locking Screw



产品名称 Product Name	材质 Material	直径(mm) Diameter	长度(mm) Length	产品编号 Product No.
梅花型中空 锁定螺钉 Torx Hollow Locking Screw	钛合金 Titanium Alloy	7.0	60、65、70、75、 80、85、90、95、 100、105、110、 115	T48870060- T48870115